




# Needlestick injuries from insulin injections are an ongoing threat to the safety of healthcare workers<sup>1</sup>



In a survey, **almost 1 in 4 nurses (~24%) experienced a needlestick injury** from giving injections to patients with diabetes<sup>1,2\*</sup>

## Consequences run far beyond the moment of accidental needlestick injury<sup>1,3,4</sup>



### Biological

Needlestick injuries put healthcare workers at risk of bloodborne pathogens, such as **HBV**, **HCV**, and **HIV**<sup>3</sup>



### Financial

Managing a needlestick injury incurs **direct and indirect costs**, with both present and future implications<sup>4</sup>



### Psychological


Needlestick injuries may leave a lasting impact on healthcare workers' **emotional well-being**<sup>1</sup>

embecta supports safe injection practices with **shorter safety insulin syringes**,<sup>†</sup> as well as **education and training resources**



# The BD SafetyGlide™ insulin syringe with a 6mm needle is designed to help protect healthcare workers from needlestick injuries

## Before injection



**Rotating safety arm** to ensure that unit scale is visible

**6mm needle** requires no pinch-up technique when given by a healthcare provider, minimizing the risk of needlestick injury through a skinfold<sup>5†</sup>

## After injection

**Sliding safety arm** enabling one-handed activation to cover the needle point after use



**Available in various barrel sizes** based on patient dosing needs (0.3mL, § 0.5mL, 1mL)

Offered in three barrel sizes:

### BD SafetyGlide™ insulin syringe

1 mL 31G × 6mm • NDC / HRI: 08290-3284-46  
• Catalog number: 328446

Quantity: 100 / Box

### BD SafetyGlide™ insulin syringe

0.5 mL 31G × 6mm • NDC / HRI: 08290-3284-47  
• Catalog number: 328447

Quantity: 100 / Box

### BD SafetyGlide™ insulin syringe

0.3 mL with ½ unit scale | 31G × 6mm  
• NDC / HRI: 08290-3284-49  
• Catalog number: 328449

Quantity: 100 / Box



## embecta is committed to being your trusted safety partner for diabetes injections

embecta offers a multi-faceted approach to needlestick injury prevention with safety-engineered devices supported by comprehensive resources, education, and training. To learn more, visit [embecta.com](https://embecta.com)

\*400 nurses in 381 different hospitals in the US voluntarily completed a survey. †Compared to >8mm. ‡Patients ≤6 years and very thin adults may require a pinch-up. §Available in half unit scale mark.

1. Frid AH, Kreugel G, Grassi G, et al. New insulin delivery recommendations. *Mayo Clin Proc.* 2016 (suppl appendix 14):1-3. 2. Lee JM, Botteman MF, Nicklasson L, Cobden D, Pashos CL. Needlestick injury in acute care nurses caring for patients with diabetes mellitus: a retrospective study. *Curr Med Res Opin.* 2015;21(5):741-747. 3. Yang L, Mullan B. Reducing the needle stick injuries in healthcare occupations: an integrative review of the literature. *ISRN Nursing.* 2011:1-11. 4. Mannocci A, De Carli G, Di Bari V, et al. How much do needlestick injuries cost? A systematic review of the economic evaluations of needlestick and sharps injuries among healthcare personnel. *Infect Control Hosp Epidemiol.* 2016; 37:635-646. 5. Frid AH, Kreugel G, Grassi G, et al. New insulin delivery recommendations. *Mayo Clin Proc.* 2016;91(9):1231-1255.

[embecta.com](https://embecta.com)

embecta, formerly part of BD. BD is the manufacturer of BD SafetyGlide insulin syringes. embecta and the embecta logo are trademarks of Embecta Corp. BD and the BD logo are trademarks of Becton, Dickinson and Company. All other trademarks are the property of their respective owners. © 2025 Embecta Corp. All rights reserved. 2505105045SGSA



Scan here for  
a demo video

